

## Platinum Electrode



KCE02, KCE03

### Platinum Coil Electrode

This electrode contains a coiled platinum wire fused in a glass tube. This electrode is used as a counter electrode in most of the electrochemical measurements.


Platinum coil electrodes are available in two different sizes and compatible with all KLYte electrochemical cell tops.

## Product details

Platinum-coil is made with a platinum wire (O.D. 0.5 mm). Part of the platinum wire is connected with copper wire inside a glass shaft to make an electrical contact. The standard ground-joint 14/23 (KCE02) and 10/18 (KCE03) enables easy assembly and are compatible with all KLYte electrochemical cell-tops.

The customized helix diameter, required for any specific application is also available.

## Electrode dimensions



KCE02, SGJ 14/23

## Specification

Electrode	KCE02	KCE03
Material	Pt	Pt
Purity	99.95%	99.95%
Shaft material	Borosilicate glass	Borosilicate glass
Electrode Plug-in-head	Compatible with KLYte banana connector cable (4mm)	Compatible with KLYte banana connector cable (2mm)
Standard Ground Joint sleeve	14/23	10/18
Shaft diameter (Top)	12mm	9mm
Shaft diameter (Bottom)	6mm	6mm
Total Length (approx.)	140mm	130mm
Pt wire Length	200mm	150mm
O.D. of Pt wire used for	0.5mm	0.5mm
Helix outer diameter	5mm (customizable)	5mm (customizable)
Height of Helix	15mm (customizable)	10mm (customizable)

## Cleaning and storing

Pure platinum metal is one of the most effective materials for counter/auxiliary electrodes due to its high electrical conductivity and resistance to corrosion. Generally, the coiled wire remains clean after the experiments, and the surface looks shiny. However, the appearance of the dull surface indicates the surface contamination. Further, surface contamination can also be detected by performing cyclic voltammetry in a pure electrolyte (e.g., 0.5mol/L aqueous  $\text{H}_2\text{SO}_4$ ). The occurrence of additional peaks other than the traditional voltammogram of platinum metal indicates the presence of surface contamination. In any case, the surface must be cleaned before using it as a counter electrode. The cleaning of the Pt-electrode can be achieved by the following methods:

➤ **The chemical method for cleaning:** Organic impurities can be cleaned with a suitable organic solvent (e.g., ethanol).

Protein deposits can be hydrolyzed with a suitable commercial enzyme-based cleaner.

Inorganic deposits can be cleaned using dilute acid and base (0.1 mol/L HCl,  $\text{HNO}_3$ , NaOH). Hot dilute acid/base solutions can be taken if the ambient temperature does not work. In general, hot 10% nitric acid removes most of the inorganic impurities.

Persistent organic/inorganic impurities on platinum surface are removed by using stronger oxidizing agents, such as freshly prepared Piranha solution (3:1 mixture of concentrated sulphuric acid and 30% hydrogen peroxide) and Aqua regia (3:1 mixture of concentrated hydrochloric acid (37%) and concentrated nitric acid). Since these oxidizing agents are strong enough to dissolve the metal from the metal surfaces, the exposure of the electrode metal into these solutions should be minimized.

*Note: One should be very careful in preparing, handling, and disposing of the Piranha solution. Mixing the solution is exothermic (hydrogen peroxide should be added slowly into concentrated sulphuric acid), and the temperature can reach above 100°C and can be explosive.*

➤ **The electrochemical method for cleaning:** Platinum electrode can be cleaned by doing multiple cyclic voltammetry in a clean solvent (10 to 20 cycles). The persistent impurities can be removed by holding the electrode either at a high oxidizing or reducing potential in dilute acid solution (0.1 mol/L sulphuric acid) for few seconds to few minutes depending upon the nature and level of impurity.

➤ **Storing:** The platinum coil of the electrode should be kept immersed in clean DI water in an airtight container while not in use.

## Optional Parts



**KEC10A**  
Banana Cable Set



**KEC10B**  
Banana Connector Pin



**KA01 (Red), KA02 (Black)**  
Alligator Clip



## Kanopy Techno Solutions

Crafting electrochemistry for lives...

### Product Information Leaflet



**Platinum Coil Electrode**  
Product ID: KCE02 and KCE03

### Contact us



[www.kanopytech.com](http://www.kanopytech.com)



[contact@kanopytech.com](mailto:contact@kanopytech.com)



+91-8604355668, +91-8004567307

A complete solution for your  
Electrochemistry research initiative...